


X8™ Genomic DNA Cartridge Kit

User Guide

FOR USE WITH:

X8™ OneTouch Instrument (Cat# X8-OT-101-IN)



Contents

Notices	3
Revision Summary	3
Introduction	4
1.1 Cartridge Kit Overview	4
1.2 Intended Use	5
1.3 Product Use Limitation	5
1.4 Storage of Eluted gDNA	5
Components & Storage	6
2.1 Cartridge Kit Components	6
2.2 Cartridge Pouch Components	6
2.3 Handling & Storage	7
X8 Workflow	9
3.1 Workflow	9
Before Each Run	10
4.1 Turn on X8 OneTouch	10
4.2 Elution Buffer Addition	10
Sample Preparation	12
5.1 User Supplied Materials	12
5.2 Cell Culture or PBMC	12
5.3 Gram-Negative Bacteria	12
5.4 Gram-Positive Bacteria	13
After Sample Preparation	14
6.1 Sample Addition	14
6.2 LyoCap Assembly	15
6.3 Starting the X8 Run	15
6.4 End of the X8 Run	15
6.5 Disposal of X8 Cartridge	15
Appendix	16
7.1 X8 Protocols	16
Troubleshooting	17

Notices

Manual Part Number

CGU-0002-EN Rev C

Legal Notices

© 2023 One BioMed Pte. Ltd. All rights reserved. Duplication and/or reproduction of all or any portion of this document without the express written consent of One BioMed Pte. Ltd., is strictly forbidden. Nothing contained herein shall constitute any warranty, express or implied, as to the performance of any products described herein. Any and all warranties applicable to any products are set forth in the applicable terms and conditions of sale accompanying the purchase of such product. One BioMed provides no warranty and hereby disclaims any and all warranties as to the use of any third-party products or protocols described herein. The use of products described herein is subject to certain restrictions as set forth in the applicable terms and conditions of sale accompanying the purchase of such product. All products and services described herein are intended FOR RESEARCH USE ONLY and NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Customer Information and Feedback

For technical information or advice, please contact our Customer Technical Support Division.

Email: support@onebiomed.com

One BioMed Pte Ltd,
4 Fusionopolis Way,
Kinesis #06
Singapore 138635
tel. +65 8189 4412

Revision Summary

Document Number

CGU-0002-EN Rev C

Effective Date

September 8, 2023

- Complete document overhaul


Introduction

1.1 Cartridge Kit Overview

The X8™ Genomic DNA Cartridge Kit is used with the X8™ OneTouch Instrument to provide a simple method for efficient, automated extraction of genomic DNA (gDNA) from cell culture, Peripheral Blood Mononuclear Cells (PBMCs) and bacterial sample types.

The ready-to-use X8 Cartridge Kits include all the necessary buffers and reagents. The X8 OneTouch Instrument is equipped with pre-programmed, automated extraction protocols, maximizing simplicity and convenience. The touchscreen Tablet PC with the intuitive and user-friendly X8 App enables users to process up to 8 samples in less than an hour. The purified gDNA can be used directly in an array of downstream applications, including Polymerase Chain Reaction, sequencing, and agarose gel electrophoresis.

Unlike typical commercial spin-column or magnetic bead technologies, the Xceler8™ Technology significantly expedites and seamlessly automates the process of gDNA extraction and purification. The novel chemistry-based approach purifies gDNA using a reversible cross-linker. It selectively binds & clusters the released gDNA from the lysed sample onto the cartridge's channel. Once impurities are washed away & sealed in the built-in waste reservoir, the purified gDNA is released and eluted using an elution buffer.



- 1.2 Intended Use** The X8 Genomic DNA Cartridge Kit (Cat# X8-GD-001-24) is intended **For Research Use Only. Not for use in diagnostic procedures.** The X8 Genomic DNA Cartridge Kit has been designed for automated extraction of gDNA from cell culture, Peripheral Blood Mononuclear Cells (PBMCs) and bacterial sample types using the X8 OneTouch Instrument (Cat# X8-OT-101-IN).
- 1.3 Product Use Limitation** The X8 Genomic DNA Cartridge Kit (Cat# X8-GD-001-24) is not intended for use with sample types other than cell culture, peripheral blood mononuclear cells (PBMCs) and bacterial samples. It is only intended for the purification of genomic DNA.
- 1.4 Storage of Eluted gDNA** If sample eluates are not processed immediately, store the eluted DNA at 4 °C.

Components & Storage

2.1 Cartridge Kit Components

The X8 Genomic DNA Cartridge Kit contains 24 Cartridges that are designed for single-use only.


The X8 Genomic DNA Cartridge Kit includes:

- 24 - X8 Genomic DNA Cartridges (PN. X8-GD01001-24)
- 5 mL - Elution Buffer (PN. X8-EB1)

2.2 Cartridge Pouch Components

The X8 Genomic DNA Cartridge Pouch includes:





- 1 - X8 Genomic DNA Cartridge (PN. X8-GD01001)
- 1 - LyoCap pouch that contains a Proteinase K and cross-linker bead




X8 Cartridge Pouch Components

2.3 Handling & Storage

In addition to the information below please follow the instructions on the X8 Cartridge Kit label for proper storage and handling of the X8 Cartridges.

SYMBOLS	DESCRIPTION
	<p>For single use only. The X8 Cartridges are intended for single use only. Do not re-use.</p>
	<p>For Research Use Only. Not for use in diagnostic procedures.</p>
	<p>Handling of infectious samples. The X8 Genomic DNA Cartridges are designed to be used with potentially infectious substances. Wear appropriate protection (e.g., gloves and goggles) when handling infectious substances. Adhere to your local and institutional guidelines for the handling and disposal of all infectious substances.</p>
	<p>Storage conditions. The X8 Genomic DNA Cartridge Kits are ethanol-free and non-volatile, therefore making the X8 Cartridges suitable for storage at ambient room temperature (15-30°C).</p>

<p>See figure below.</p>	<p>Storage and handling orientations. Improper shipping, storage or handling orientations may result in reagent blockage within the X8 Cartridges, which can then adversely impact the extraction performance.</p> <p>Upon receiving the X8 Cartridges-</p> <ul style="list-style-type: none">• Rest the X8 Cartridges in their proper storage orientation overnight prior to use. <p>Before commencing the nucleic acid extraction-</p> <ul style="list-style-type: none">• Rest the X8 Cartridges in its proper storage orientation for a minimum of 30 min prior to use.
--------------------------	--




Optimal storage and handling orientation.

Set up on a flat surface with the label-side facing up.

For more information on safe handling, please refer to the Product Safety Data Sheets (SDS, CGU-0003) available at www.onebiomed.com/products.

X8 Workflow

3.1 Workflow




Before Each Run

4.1 Turn on X8 OneTouch


Ensure the X8 OneTouch Instrument is turned on and equilibrated. Please refer to the X8 OneTouch Operating Manual (CGU-0001) for additional details.

4.2 Elution Buffer Addition

1. Orient the X8 Cartridge on a flat surface with the reagent wells facing the user.




2. Gently lift the blue sticker tab on the X8 Cartridge and pierce the foil at the marked position using a sterile pipette tip.



- Carefully aspirate 150 μ L of the DNA Elution Buffer (PN. X8-EB1), ensuring there are no air bubbles.

Note. Please refer to the appendix for additional details on varying volumes of the Elution Buffer.

- Gently insert the pipette tip into the pierced foil until it touches the base of the well. Then position the pipette tip at a 45° angle to touch the front base of the well and then slowly dispense the DNA Elution Buffer into the well without introducing any air bubbles.




Note. Incorrect addition of the Elution Buffer into the X8 Cartridge can result in little or no elution of the purified nucleic acid.

- Peel the backing paper on the tape and apply the tape to seal the pierced foil. Run your finger across the tape to ensure that it is tightly sealed.

Note. Verify the foil is clean and dry before applying the tape. Ensure that there is no visible air gap between the pierced foil and the tape. Improper application of the tape can result in sample cap error.

Note. Do not shake or tilt the X8 cartridge once the Elution Buffer has been added.



Sample Preparation

5.1 User Supplied Materials

1. Micropipettes and sterile pipette tips.
2. Sterile microcentrifuge tubes (0.5 mL or 1.5 mL).
3. 1X PBS.
4. Centrifuge (with adapter for 1.5 mL / 2 mL microcentrifuge tubes).

For Gram-positive bacteria extraction (if applicable):

5. Lysozyme reconstituted as a 25 mg/mL stock in lysozyme buffer (20 mM Tris-HCl pH 8.0, 2 mM EDTA, 1.2% TritonX-100).
6. Heating block.

5.2 Cell Culture or PBMC

1. Obtain the cells of interest in suspension.
2. Centrifuge to pellet the cells and discard supernatant. Recommended speed: 200 x g for 5 minutes.
3. Wash the cells by adding at least 1 mL of PBS per 1×10^6 of cells and mix gently by pipetting to resuspend cells. Repeat step 2.
4. Resuspend the cell pellet to the final concentration of $\sim 1 \times 10^6$ cells in 200 μ L of PBS.

Note. For PBMC, resuspend the cell pellet to the final concentration of 2×10^6 cells in 200 μ L of PBS.

5.3 Gram-Negative Bacteria

1. Obtain bacteria from an overnight culture.
2. Centrifuge to pellet the bacteria and discard supernatant. Recommended speed: 6080 x g for 5 minutes.
3. Wash the bacteria by adding at least 1 mL of PBS per 1×10^9 of bacteria and mix gently by pipetting to resuspend bacteria. Repeat step 2.
4. Resuspend the bacteria pellet to the final concentration of $\sim 1 \times 10^9$ bacteria in 200 μ L of PBS.

5.4 Gram-Positive Bacteria

1. Obtain bacteria from an overnight culture.
2. Centrifuge to pellet the bacteria and discard supernatant. Recommended speed: 6080 x g for 5 minutes.
3. Wash the bacteria by adding at least 1 mL of PBS per 1×10^9 of bacteria and mix gently by pipetting to resuspend bacteria. Repeat step 2.
4. Resuspend the bacteria pellet to the final concentration of $\sim 1 \times 10^9$ bacteria in 80 μ L of PBS.
5. Add 120 μ L of 25 mg/mL lysozyme in 20 mM Tris-HCl (pH 8.0), 2 mM EDTA, 1.2% Triton X-100.
6. Incubate at 37 °C for 30 minutes on a heat block.


Note. Proceed to section 6.1 immediately after lysozyme treatment.

After Sample Preparation

6.1 Sample Addition


1. Open the sample well cap and transfer all of the prepared sample volume (as per section 5) into the sample well.

Note. Ensure the pipette tip does not touch the rim of the sample well or the air hole within the sample well.



2. Once the sample has been added, close the sample well cap tightly by pushing the cap down for at least 1 second.


Note. In case of a "Check sample cap" error message, take out the X8 Cartridge and re-check to make sure the sample cap is sealed tightly before inserting the Cartridge into the bay.



6.2 LyoCap Assembly

1. Next open the LyoCap pouch and retrieve the LyoCap. Press the LyoCap (foil side down) onto the capsule holder on the X8 Cartridge.

Note. There is no specific orientation to install the LyoCap. Minimize the time between LyoCap assembly and starting X8 run.



6.3 Starting the X8 Run

1. The X8 Cartridge is now ready to use with the X8 OneTouch Instrument.
2. Please place a sterile microcentrifuge tube in the tube holder before inserting the X8 Cartridge.
3. For the operation of X8 OneTouch Instrument, please refer to the X8 OneTouch Operating Manual (CGU-0001).
4. For the selection of protocols based on the sample types, please refer to the Appendix for the X8 Protocols.

6.4 End of the X8 Run

1. When the run is completed, the LED indicator will turn BLUE and a "Done" message will show on the X8 App.
2. Remove the X8 Cartridge from the bay.
3. Carefully remove the microcentrifuge tube from the tube holder and close the cap tightly.

Note. If sample eluates are not processed immediately, store the eluted DNA at 4 °C.

6.5 Disposal of X8 Cartridge

1. Treat the used X8 Cartridge as biohazard and discard as per local and organizational guidelines. Please refer to X8 Genomic DNA Cartridge Kit Safety Data Sheet (CGU-0003).

Appendix

7.1 X8 Protocols

SAMPLE TYPE	TOTAL CELL CONC. INPUT	PROTOCOL NAME	ELUTION VOLUME
Cultured cells	1×10^6	GD_CC_V#_#	150 - 200 μ L
Gram Negative	1×10^9	GD_GN_V#_#	150 - 200 μ L
<i>E.coli</i>	1×10^9	GD_EC_V#_#	150 - 200 μ L
Gram Positive	1×10^9	GD_GP_V#_#	150 - 200 μ L
Cultured cells	$> 0.1 \times 10^6$ to $< 1 \times 10^6$	GD_CC_V#_#_ Low_EB	50 - 100 μ L
Gram Negative	$> 0.1 \times 10^9$ to $< 1 \times 10^9$	GD_GN_V#_#_ Low_EB	50 - 100 μ L
<i>E.coli</i>	$> 0.1 \times 10^9$ to $< 1 \times 10^9$	GD_EC_V#_#_ Low_EB	50 - 100 μ L
Gram Positive	$> 0.1 \times 10^9$ to $< 1 \times 10^9$	GD_GP_V#_#_ Low_EB	50 - 100 μ L

1. For cultured cells sample concentration between 0.1×10^6 to less than 1×10^6 , or bacterial cells between 0.1×10^9 to less than 1×10^9 , please use Low EB protocol and Elution Buffer volume of 50 to 100 μ L.
2. The "...V#_#" in the protocol name refers to the version of the protocol.

Troubleshooting

For additional technical information or advice, please contact our Customer Technical Support Division.

Email: support@onebiomed.com

TECHNICAL PROBLEM	DESCRIPTION
<p>“Check sample cap” error</p>	<p>Reapply the sample cap firmly.</p> <p>Ensure that there is no visible air gap between the pierced foil and the tape.</p>
<p>Lower than expected DNA yield</p>	
<p>Incorrect sample concentration</p> <p>Incorrect X8 protocol used</p>	<p>The genomic DNA yield from samples depends on the sample concentration and the corresponding X8 protocol.</p> <p>Refer to Appendix for the recommended sample concentration and X8 protocol.</p>
<p>DNA degradation after storage</p>	<p>It is recommended to store the eluted DNA at 4 °C.</p>
<p>Sample did not fully enter the channel</p>	<p>Maximum 200 µL of sample volume is to be added, ensure there are no air bubbles while loading into the sample well.</p> <p>Avoid touching the airhole while loading into the sample well.</p>

Lower than expected purity ratios	
Sample overload	Refer to Appendix for the recommended sample concentration and X8 protocol.
Compromised washing	Improper storage or handling of the X8 Cartridges may result in reagent blockage in the X8 Cartridge.
Viscous eluate	Ensure that no more than 2×10^6 cultured cells or 1×10^9 bacterial cells are used per sample.
Elution Buffer related issues	
Too much elution	Improper storage or handling orientations may result in reagent blockage in the X8 Cartridge.
Too little or no elution	<p>Incorrect addition of the elution buffer. Please refer to section 4.2 for Elution Buffer Addition.</p> <p>There will be dead volume retained in the X8 Cartridge channel after elution of nucleic acid.</p>

RESOURCES

We are dedicated to helping you get the most out of your Xceler8 Platform by offering multiple helpful resources:

Products page

Explore our ever-growing portfolio of products that Xceler8 your sample preparation workflow.

onebiomed.com/products/

Support

Email us for product support and feedback.

support@onebiomed.com

Contact us

[web onebiomed.com](https://onebiomed.com) | [email info@onebiomed.com](mailto:info@onebiomed.com)

address 4 Fusionopolis Way, Kinesis, #06, Singapore 138635

Copyright © One BioMed Pte. Ltd. All rights reserved.
FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.
CGU-0002-EN Rev C

The logo for One BioMed, featuring the word "One" in a dark blue font and "BioMed" in an orange font, with a stylized orange and blue circular graphic element above the "i" in "BioMed".

One BioMed